

Product datasheet for **MR223661L4V**

Atp2b2 (NM_001036684) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Atp2b2 (NM_001036684) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Atp2b2
Synonyms:	D6Abb2e; dfw; Gena300; jog; PMCA2; Tmy; wms; wri
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001036684
ORF Size:	3594 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223661).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001036684.2
RefSeq Size:	6925 bp
RefSeq ORF:	3597 bp
Locus ID:	11941
UniProt ID:	Q9R0K7
Cytogenetics:	6 52.85 cM
Gene Summary:	This magnesium-dependent enzyme catalyzes the hydrolysis of ATP coupled with the transport of calcium out of the cell. Plays a role in maintaining balance and hearing. [UniProtKB/Swiss-Prot Function]


[View online »](#)