

Product datasheet for **MR223683L4V**

Atp1b3 (NM_007502) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Atp1b3 (NM_007502) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Atp1b3
Synonyms:	AA409958; AI664000; AW212096
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007502
ORF Size:	834 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223683).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007502.4
RefSeq Size:	1833 bp
RefSeq ORF:	837 bp
Locus ID:	11933
UniProt ID:	P97370
Cytogenetics:	9 50.31 cM
Gene Summary:	This is the non-catalytic component of the active enzyme, which catalyzes the hydrolysis of ATP coupled with the exchange of Na(+) and K(+) ions across the plasma membrane. The exact function of the beta-3 subunit is not known.[UniProtKB/Swiss-Prot Function]



[View online »](#)