

Product datasheet for MR223723L3V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Abcc8 (NM_011510) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: Abcc8 (NM_011510) Mouse Tagged ORF Clone Lentiviral Particle

Symbol: Abcc8

Synonyms: D930031B21Rik; Sur; SUR1

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

 Tag:
 Myc-DDK

 ACCN:
 NM_011510

 ORF Size:
 4764 bp

ORF Nucleotide

Sequence:

The ORF insert of this clone is exactly the same as(MR223723).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through

naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: <u>NM 011510.3</u>, <u>NP 035640.2</u>

 RefSeq Size:
 4952 bp

 RefSeq ORF:
 4767 bp

 Locus ID:
 20927

Cytogenetics: 7 29.66 cM

