

## Product datasheet for **MR223733**

### **Il1rap (NM\_001159317) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Il1rap (NM_001159317) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Il1rap
Synonyms:	6430709H04Rik; AI255955; AV239853; IL-1RAcP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MR223733 representing NM\_001159317  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGGACTTCTGTGGTATTTGATGAGTCTGTCTTCTATGGGATCCTGCAGAGTCATGCTTCGGAGCGCT  
GTGATGACTGGGACTAGATACCATGCGACAAATCCAAGTGTGTTGAAGATGAGCCGGCTCGAATCAAGTG  
CCCCCTCTTTGAACACTTCTGAAAGTACAACACTACAGCACTGCCATTCTCTGGCCTTACCTGATCTGG  
TACTGGACCAGGCAAGACCGGGACCTGGAGGAGCCATTAACCTCCGCCTCCAGAGAATCGCATCAGTA  
AGGAGAAAAGATGTGCTCTGGTCCGGCCACCCTCTCAATGACACGGGCAATTACACCTGCATGTTGAG  
GAACACAACCTTACTGCAGCAAAGTTGCATTTCCCTGGAAGTTGTTCAGAAGGACAGCTGTTTCAATTCT  
GCCATGAGATCCAGTGCACAAGATGTATATTGAACATGGCATTATAAGATCACATGTCCAAATGTAG  
ACGGATACTTTCCCTCCAGTGTCAAACCATCGTCACTTGGTATAAGGGTGTACTGAAATAGTGGACTT  
TCATAATGTACTACCCGAGGGCATGAACTTGAGCTTTTTTCATCCCCTTGGTTTCAAATAACGGAAATTAC  
ACATGTGTGGTTACATATCCTGAAAACGGACGTCTCTTTCACCTACCAGGACTGTGACTGTAAAGGTGG  
TGGGCTACCAAAGGATGCATTGCCACCCAGATCTATTCTCAAATGACCGTGTGTCTATGAGAAAGA  
ACCAGGAGAGGAACTGGTTATCCCTGCAAAGTCTATTTCAAGTTTCATTATGGACTCCCAATGAGGTC  
TGGTGGACCATTGATGAAAAGAAGCCTGATGACGTACAGTCGACATCACTATTAATGAAAGTGAAGTT  
ATCTTCAACGGAAGATGAAACAAGGACTCAGATTTTGAGCATCAAGAAAGTCAACCCGGAGGATCTCAG  
GCGCAACTATGTCTGTATGCTCGAAATACCAAAGGGGAAAGTGAAGCAGGCTGCCAAGGTGAAACAGAAA  
GTCATACCACCAAGGTACACAGTAGAACTCGCCTGTGGTTTTGGAGCCACGGTCTTCTGGTAGTGGTTC  
TCATTGTGGTTTTACCATGTTTACTGGCTGGAGATGGTCTCTTTTACCGAGCTCACTTTGGAACAGATGA  
AACAAATCTTGATGAAAAGGAGTATGATATTTATGTTTCTATGCAAGAAATGTGGAAGAAGAGGAATTT  
GTGCTGCTGACGCTGCGTGGAGTTTTGGAGAATGAGTTTGGATACAAGCTGTGCATCTTCGACAGAGACA  
GCCTGCCTGGGGGAAATACCGTGGAAAGCAGTTTTTGTATTCATTAGAGGAGCCGAAGGATGATTGTTGT  
CCTGAGCCCTGACTATGTGACAGAAAAGAGCATCAGCATGCTGGAGTTTAAAGCTGGGTGTCATGTGCCAG  
AACTCCATTGCCACTAAGCTCATTGTGGTGGAGTACCGTCCGCTTGAGCAACCCCATCCAGGCATCATGC  
AGCTGAAGGAGTCTGTGCTTTTGTAAAGCTGGAAGGGAGAAAAGTCCAAACATTCTGGCTCCAAGTTCTG  
GAAGGCCTTGCGTTTGGCTCTTCCCTGAGAAGTCTGAGCGCCAGCTCCGGCTGGAATGAGAGCTGTTCT  
TCTCAGTCTGACATCAGTCTGGATCATGTTGAGAGGAGAAGTCGTTTGAAGAGCCCCAGAACTCCGAA  
GCTCAGAGAGGGTGTCTGGAGCAGAGCCAGCCCCGGGCACGATGTCCAAGCACCGAGGAAACCCTCAGC  
AGCCTGTGCTGTGTGTACCTACTGTGAAGGAGAAAAGTCACTCAGGAGCAAGAGCCGGGCAGAGATG  
CACACGCATCCCCAGTGGGAAACACACCTC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR223733 representing NM\_001159317  
Red=Cloning site Green=Tags(s)


MGLLWYLSLSFYGILQSHASERCDDWGLDTMRQIQVFEDEPARIKCPLFEHFLKYNSTAHSSGLTLIW  
 YWTRQDRDLEEFINFRLENRISKEKDVLFWRPTLLNDTGNYTCLRNNTYCSKVAFPLEVVQKDSFCNS  
 AMRFPVHKMYIEHGIHKITCPNVDGYFPSSVKPSVTWYKGCIEIVDFHNVLP EGMNLSFFIPLVSNNGNY  
 TCVVYTPENGRFLHLTRVTVKVVGSPKDALPPQIYSPNDRVVYEKEPGEELVIPCKVYFSFIMDSHNEV  
 WWTIDGKKPDDVTVDITINESVSYSTEDETRITQILSIKKVTPEDLRNNYVCHARNTKGEAEQAQAKVKQK  
 VIPPRYTVELACGFGATVFLVVVLI VVYHVVWLEMLVLFYRAHFGTDEITLDGKEYDIYVSYARNVEEEEF  
 VLLTLRGVLENEFGYKLCIFDRDSLPGGNTVEAVDFDIQRSRMIVVLSPDYVTEKSISMLEFKLGVMCQ  
 NSIATKLI VVEYRPLEQPHPGIMQLKESVSFVSWKGEKSKHSGSKFWKALRLALPLRSLASSGWNEESC  
 SQSDISLDHVQRRSRLKEPPELRSSERVSGAEPAPGTMSKHRGKPSAACRCCVTYCEGESHLRSKRAEM  
 HTHPQWETHL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9005\\_a12.zip](https://cdn.origene.com/chromatograms/mm9005_a12.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001159317

**ORF Size:** 1920 bp


**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001159317.1</a></u> , <u><a href="#">NP_001152789.1</a></u>
<b>RefSeq Size:</b>	4031 bp
<b>RefSeq ORF:</b>	1923 bp
<b>Locus ID:</b>	16180
<b>Cytogenetics:</b>	16 B2
<b>MW:</b>	74 kDa
<b>Gene Summary:</b>	<p>Coreceptor for IL1RL2 in the IL-36 signaling system. Coreceptor with IL1R1 in the IL-1 signaling system. Associates with IL1R1 bound to IL1B to form the high affinity interleukin-1 receptor complex which mediates interleukin-1-dependent activation of NF-kappa-B and other pathways. Signaling involves the recruitment of adapter molecules such as TOLLIP, MYD88, and IRAK1 or IRAK2 via the respective TIR domains of the receptor/coreceptor subunits. Recruits TOLLIP to the signaling complex. Does not bind to interleukin-1 alone; binding of IL1RN to IL1R1, prevents its association with IL1R1 to form a signaling complex. The cellular response is modulated through a non-signaling association with the membrane IL1R2 decoy receptor. Secreted forms (isoforms 2 and 3) associate with secreted ligand-bound IL1R2 and increase the affinity of secreted IL1R2 for IL1B; this complex formation may be the dominant mechanism for neutralization of IL1B by secreted/soluble receptors. Coreceptor for IL1RL1 in the IL-33 signaling system. Can bidirectionally induce pre- and postsynaptic differentiation of neurons by trans-synaptically binding to PTPRD (PubMed:25908590). May play a role in IL1B-mediated costimulation of IFNG production from T-helper 1 (Th1) cells (By similarity). [UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for MR223733