

## Product datasheet for **MR223734L4V**

### Adgre5 (NM\_011925) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Adgre5 (NM_011925) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Adgre5
Synonyms:	AA409984; Cd97; TM7LN1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011925
ORF Size:	2454 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223734).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011925.2</a> , <a href="#">NP_036055.2</a>
RefSeq Size:	3253 bp
RefSeq ORF:	2457 bp
Locus ID:	26364
Cytogenetics:	8 40.22 cM
Gene Summary:	Receptor potentially involved in both adhesion and signaling processes early after leukocyte activation. Plays an essential role in leukocyte migration.[UniProtKB/Swiss-Prot Function]



[View online »](#)