

Product datasheet for MR223737

Dusp22 (NM_134068) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dusp22 (NM_134068) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Dusp22
Synonyms: JKAP; JSP1; LMW-DSP2; MKPX
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR223737 representing NM_134068
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGGGAGTGGGATGAGCCAGATCCTGCCGGCCCTGTACATTGGCAACTCAAAGACGCAAGAGATGCAG
 AACAGTTGAGCAGGAACAAGGTGACACACATTCTTTCTGTGCACGATACTGCCAGGCCCATGTTGGAGGG
 AGTTAAATACCTGTGTATCCAGCGGCAGACACACCATCTCAAAACCTGACAAGACATTTCAAAGAAAGC
 ATTAATTCATTCATGAGTGCCGACTCCAGGGTGAGAGCTGTCTTGTACATTGCCTGGCTGGGTCTCCA
 GGAGTGTGACATTGGTGATCGCATACATCATGACTGTCACCGACTTTGGCTGGGAAGATGCCTTGCACAC
 TGTTCTGTCGGGGAGGTCTGTGCCAACCCCAACTGGGCTTTCAAAGGCAGCTGCAGGAGTTTGAGAAA
 CATGAAGTGCACCAAGTATCGGCAATGGCTGAGAGAAGAGTATGGAGAGAACCCTTTGCGGGATGCAGAAG
 AAGCAAAAATATTCTGGCTGCCCGGGAATTCGAAGTACTGGGCTTTCTCAGAAGACTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR223737 representing NM_134068
 Red=Cloning site Green=Tags(s)

MGSGMSQILPGLYIGNFKDARDAEQLSRNKVTHILSVHDTARPMLEGVKYL CIPAAADTPSQLTRHFKEK
 IKFIEHCR LQGESCLVHCLAGVSRVTLVIAYIMTVDFGWEDALHTVRAGRSCANPNLGFQRQLQEFEK
 HEVHQYRQWLREEYGENPLRDAEEAKNILAAPGILKYWAF LRL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



Cloning Scheme:



ACCN: NM_134068

ORF Size: 552 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_134068.3](#), [NP_598829.1](#)

RefSeq Size: 1058 bp

RefSeq ORF: 555 bp

Locus ID: 105352


UniProt ID: [Q99N11](#)

Cytogenetics: 13 A3.2

MW: 21.4 kDa

Gene Summary: Activates the Jnk signaling pathway. Dephosphorylates and deactivates p38 and stress-activated protein kinase/c-Jun N-terminal kinase (SAPK/JNK).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR223737