

Product datasheet for **MR223740L4V**

Abra (NM_175456) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Abra (NM_175456) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Abra
Synonyms:	C130068O12Rik; STARS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_175456
ORF Size:	1125 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223740).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_175456.4 , NP_780665.1
RefSeq Size:	1762 bp
RefSeq ORF:	1128 bp
Locus ID:	223513
UniProt ID:	Q8BUZ1
Cytogenetics:	15 B3.1
Gene Summary:	Acts as an activator of serum response factor (SRF)-dependent transcription possibly by inducing nuclear translocation of MKL1 or MKL2 and through a mechanism requiring Rho-actin signaling.[UniProtKB/Swiss-Prot Function]



[View online »](#)