

Product datasheet for MR223754L3

Mag (NM_010758) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mag (NM_010758) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Mag
Synonyms:	Gm; Gma; sigle; siglec-4a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223754).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	


Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Plasmid Map:


ACCN: NM_010758

ORF Size: 1749 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_010758.2](#), [NP_034888.1](#)

RefSeq Size: 2427 bp

RefSeq ORF: 1749 bp

Locus ID: 17136

UniProt ID: [P20917](#)

Cytogenetics: 7 19.26 cM

Gene Summary: This gene encodes a type I membrane protein and member of the immunoglobulin-like superfamily. It is expressed in myelinating glial cells, including oligodendrocytes of the central nervous system and Schwann cells of the peripheral nervous system. Mice lacking the encoded protein express abundant myelin, but suffer long-term axon degeneration, altered distribution of channels and adhesion molecules at nodes of Ranvier, and altered axon cytoskeletal structure. While not required for myelination, the encoded protein enhances axon-myelin stability, helps to structure nodes of Ranvier, and regulates the axon cytoskeleton. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2016]