

## Product datasheet for **MR223755L3V**

### Cep76 (NM\_001081073) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cep76 (NM_001081073) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cep76
Synonyms:	6230425F05Rik; 9630013E05
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001081073
ORF Size:	1980 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223755).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001081073.1</a> , <a href="#">NP_001074542.1</a>
RefSeq Size:	3979 bp
RefSeq ORF:	1980 bp
Locus ID:	225659
UniProt ID:	<a href="#">Q0VEJ0</a>
Cytogenetics:	18 E1
Gene Summary:	Centrosomal protein involved in regulation of centriole duplication. Required to limit centriole duplication to once per cell cycle by preventing centriole reduplication (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)