

Product datasheet for **MR223765L3V**

Atp10d (NM_153389) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Atp10d (NM_153389) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Atp10d |
| Synonyms: | 9830145H18Rik; ATPVD; D5Buc24e |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_153389 |
| ORF Size: | 4407 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR223765). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_153389.3 , NP_700438.3 |
| RefSeq Size: | 6021 bp |
| RefSeq ORF: | 4410 bp |
| Locus ID: | 231287 |
| UniProt ID: | Q8K2X1 |
| Cytogenetics: | 5 C3.2 |



[View online »](#)

Gene Summary:

Catalytic component of a P4-ATPase flippase complex which catalyzes the hydrolysis of ATP coupled to the transport of aminophospholipids from the outer to the inner leaflet of various membranes and ensures the maintenance of asymmetric distribution of phospholipids. Phospholipid translocation seems also to be implicated in vesicle formation and in uptake of lipid signaling molecules (By similarity).[UniProtKB/Swiss-Prot Function]