

Product datasheet for **MR223795L4V**

Tslp (NM_021367) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tslp (NM_021367) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tslp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021367
ORF Size:	420 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223795).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021367.2
RefSeq Size:	1143 bp
RefSeq ORF:	423 bp
Locus ID:	53603
UniProt ID:	Q9JIE6
Cytogenetics:	18 18.12 cM
Gene Summary:	Cytokine that induces the release of T-cell-attracting chemokines from monocytes and, in particular, enhances the maturation of CD11c(+) dendritic cells. Can induce allergic inflammation by directly activating mast cells (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)