

Product datasheet for MR223821L1

Gm7325 (Mymx) (NM_001177469) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gm7325 (Mymx) (NM_001177469) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Mymx
Synonyms:	EG653016; Esgp; Gm7325; minion; myomerge
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223821).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_001177469
ORF Size:	252 bp




[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001177469.1 , NP_001170940.1
RefSeq Size:	811 bp
RefSeq ORF:	255 bp
Locus ID:	653016
UniProt ID:	Q2Q5T5
Cytogenetics:	17 B3
Gene Summary:	Myoblast-specific protein that mediates myoblast fusion, an essential step for the formation of multi-nucleated muscle fibers (PubMed:28386024, PubMed:28569745, PubMed:28569755, PubMed:30197239). Involved in membrane fusion downstream of the lipid mixing step mediated by MYMK (PubMed:30197239). Acts by generating membrane stresses via its extracellular C-terminus, leading to drive fusion pore formation (PubMed:30197239). Acts independently of MYMK (PubMed:30197239). Involved in skeletal muscle regeneration in response to injury by mediating the fusion of satellite cells, a population of muscle stem cells, with injured myofibers (PubMed:29581287).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR223821L1



Double digestion of MR223821L1 using Sgf1 and MluI