

Product datasheet for **MR223831**

Myt1l (NM_001093776) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Myt1l (NM_001093776) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Myt1l
Synonyms:	2900046C06Rik; 2900093J19Rik; C630034G21Rik; mKIAA1106; Nztf1; Pmng1; Png-1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>MR223831 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)


TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACGTGGACTCTGAGGAGAAGCGCCATCGCACACGGTCCAAAGGGGTTCGAGTTCCTGTGGAGCCAG
CCATACAAGAGCTGTTCCAGCTGTCCACTCCAGGCTGCGACGGCAGTGGTCACGTCAGTGGCAAATATGC
ACGACACAGAAGTGTATATGGTTGTCCCTTGGCTAAAAAAGAAAAACGCAAGATAAACAGCCCCAAGAA
CCTGCTCCAAGCGAAAACCATTTGCAGTAAAAGCAGATAGTTCCCTCAGTAGACGAATGTTATGAGAGTG
ATGGTACTGAAGACATGGATGATAAGGAGGAAGATGATGATGAGGAGTTCCTGAAGACAATGATGAGCA
AGGGGATGATGACGACGAAGATGAGGTGGATCGGGAAGACGAGGAGGAGATCGAGGAGGAAGATGATGAA
GAAGATGATGATGATGAAGATGGTGACGATGTAGAAGAGGAAGAAGAGGATGATGATGAAGAGGAGGAAG
AAGAGGAAGAGGAAGAAGAAAAATGAAGACCATCAAATGAGTTGTACTCGAATAATGCAGGACACAGACAA
GGATGATAACAACAATGATGAGTATGATAACTATGATGAACTGGTAGCTAAGTCGCTATTAATCTTGGC
AAAATTGCTGAGGATGCAGCATACCGAGCCAGGACTGAATCAGAGATGAACAGCAATACCTCAAATAGTC
TGGAGGACGATAGTGACAAAAACGAAAACCTCGGTGCGAAAAGCGAACTGAGTCTAGACTAGACAGTGA
TGTTGTTAGAGAAACAGTGGACTCCCTTAAGCTGTTAGCAACAAGGACATGGTGTGTGCTATCAGAGAAT
ATCAGTGACAGAAGTTATGCTGAGGGGATGTCACAGCAGGACAGTAGAAAATATGAACTATGTCATGCTAG
GGAAGCCCATGAACAATGGACTCATGGAGAAGATGGTGGAGGAGAGTGTGAGGAAGTGTGTCTAAGTAG
TCTAGAGTGCCTGAGGAACCAAGTCTTTGACCTGGCCAGGAACTCAGCGAGACCAACCCACAGGACAGG
AGTCAGCCACCCAACATGAGTGTGCGCCAACATGTCCGGCAAGAGGACGACTTCCCTGGGAGGACGCCAG
ACAGGAGTACTCGGATATGATGAACCTTATGCGGCTGGAGGAGCAGCTCAGTCCAGGTCTAGAACGTT
CTCCAGCTGTGCAAGGAGGATGGGTGTCATGAGAGGGATGATGACACCCTCAGTGAACCTCAGACAGG
TCTGAGGAAGTGTGACATGACCAAGGGCAACCTGACTCTGCTAGAGAAAGCCATTGCCTTGGAGACAG
AGAGAGCCAAGGCCATGCGGGAGAAGATGGCCATGGATGCTGGGAGAAGGATAACCTGAGATCCTATGA
GGACCAGTCTCCAAGACAGCTGGCTGGGGAAGACAGAAAATCAAATCCAGTGACAGCCATGTCAAAAAG
CCATACTATGATCCCTCAAGAACAGAAAAGAGAGAGAGCAAGTGTCCAACCCCGGGTGTGATGGAACCG



GCCACGTAACGGGCTTTACCCGCATCACCGCAGTCTGTCTGGATGCCCGCACAAAGATAGGGTCCCTCC
 AGAAATCTTGCCATGCATGAAAATGTTCTCAAGTGTCCCACTCCAGGCTGCACAGGGCGAGGGCATGTG
 AATAGCAACAGGAACTCGCACAGAAGCCTCTCTGGATGCCCAATTGCTGCTGCAGAAAACTGGCAAAGG
 CCCAAGAGAAACACCAGAGCTGTGATGTGTCCAAATCCAACCAGGCTCAGACCGAGTCTCAGGCCAAT
 GTGCTTTGTCAAACAGCTTGAGATTCTCAGTATGGCTACAGAAACAATGTTCCCAACACACACCGC
 TCCAACCTGGCCAAGGAGCTTGAGAAATACTCCAAGACTTCGTTTGAGTACAACAGTTACGACAACCATA
 CTTATGGCAAAAAGAGCCATAGCTCCCAAGGTGCAAACCAGGGACATATCCCCAAAGGATATGACGATGC
 CAAGCGGTACTGCAAGAATGCCAGCCCCAGCAGCAGCACCACCAGCAGCTATGCACCTAGCAGCAGCAGC
 AACCTCAGCTGTGGTGGTGGCAGCAGCGCCAGTAGCACGTGTAGCAAGAGCAGCTTTGACTACACACATG
 ACATGGAGGCCGCACACATGGCAGCCACAGCCATTCTCAACCTGTCCACACGTTGTGCGTAAATGCCACA
 GAACCTGTCCACCAAGCCACAGGACCTGTGACTGCCCGAAACCCAGACATGGAGGTGGATGAGAATGGC
 ACCCTGGACCTGAGCATGAACAAGCAGAGGCTCGAGACAGCTGTGCCAGTCTGACACCCCTGGAAC
 CCATGTCTCCGCAGCAGGCGGTGATGAGCAGCCGATGCTTCCAGCTGAGCGAGGGGGATTGCTGGGA
 CTTGCCTGTAGACTACACAAAAATGAAGCCTCGGAGGGTAGATGAGGATGAGCCAAAGAGATTACCCCA
 GAAGACTTGGACCAATCCAGGAGGCTCTGGAAGAAAGACGGTATCCAGGGGAGGTGACCATCCCAAGCC
 CCAAAACCAAGTACCCTCAGTGCAAGGAAAGCAAAAAGGACTTAATAACTCTGTCTGGCTGCCCCCTGGC
 GGACAAAAGCATTGCAAGTATGCTGGCCACCAGTCCCAAGAGCTCAAGTGCCCCACCCCTGGCTGTGAC
 GGTCTGGACACATCACTGGCAATTACGCTTCTCATCGAAGCCTTCTGGGTGCCGAGAGCAAAGAAGA
 GTGGCATCCGGATAGCACAGAGCAAAGAGGACAAGGAAGACCAGGAGCCAATCAGGTGTCCGGTACCTGG
 CTGTGACGGTCAGGGACACATCACTGGGAAGTATGCATCCACCAGCAGCGCTCCGGGTGTCCCTTGGCA
 GCCAAGAGGCAGAAAGATGGGTACCTTAATGGCTCCCAGTCTCCTGGAAGTCGGTCAAGACGGAGGGCA
 TGTCTGCCCTACCCCGGGTGTGATGGGTGAGGACACGTGAGTGGCAGCTTCTCACACCCGAGCTT
 GTCAGGATGTCCAAGAGCCACATCAGCAATGAAGAAAGCAAAGCTGTCTGGAGAACAGATGTTGACTATC
 AAGCAGCGAGCCAGCAACGGTATAGAAAATGATGAAGAAATCAAGCAGTTAGATGAAGAGATCAAGGAGC
 TTAATGAGTCCAATCCAGATGGAGGCTGACATGATCAAACCTCAGAACTCAGATCACCACAATGGAGAG
 CAACCTGAAGACGATTGAGGAGGAGAACAAGTCAATGAACAGCAGAATGAGTCGCTCTTGACAGGTTG
 GGCAACCTGAGCCAGTCCCTGATCCACAGCCTCGCCAACATCCAGCTGCCTCACATGGATCCAATCAATG
 AACAAAATTTTGATGCTTACGTGACTACTTTGACGGAAATGTATACAAATCAAGATCGTTATCAGAGTCC
 AGAAAATAAGCCCTACTGGAAAATATAAAGCAGGCTGTGAGAGGAATTCAGGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Plasmid Map:


ACCN: NM_001093776

ORF Size: 3555 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001093776.1](#), [NP_001087245.1](#)

RefSeq Size: 7104 bp

RefSeq ORF: 3558 bp

Locus ID: 17933

Cytogenetics: 12 11.86 cM

MW: 132.8 kDa

Gene Summary:

Transcription factor that plays a key role in neuronal differentiation by specifically repressing expression of non-neuronal genes during neuron differentiation (PubMed:28379941). In contrast to other transcription repressors that inhibit specific lineages, mediates repression of multiple differentiation programs (PubMed:28379941). Also represses expression of negative regulators of neurogenesis, such as members of the Notch signaling pathway, including HES1 (PubMed:28379941). The combination of three transcription factors, ASCL1, POU3F2/BRN2 and MYT1L, is sufficient to reprogram fibroblasts and other somatic cells into induced neuronal (iN) cells in vitro (PubMed:20107439, PubMed:24243019, PubMed:27281220). Directly binds the 5'-AAGTT-3' core motif present on the promoter of target genes and represses transcription by recruiting a multiprotein complex containing SIN3B (PubMed:28379941). The 5'-AAGTT-3' core motif is absent from the promoter of neural genes (PubMed:28379941).[UniProtKB/Swiss-Prot Function]