

## Product datasheet for MR223833

### Pygo1 (NM\_028116) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Pygo1 (NM_028116) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pygo1
Synonyms:	2600014C22Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR223833 representing NM_028116 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCAGCGGAACAGGACAAGGAGCCCATCGCGCTGAAGAGAGTTAGAGGTGGTGACAGTGGACTGGATG  
GGTTAGGAGGGCCCAATATACAACCTCGGAAGCCAGATAAGAAAAACGCAAGGCCAACACACAGGGATC  
TTCCTTTCCTCCGCTGTCAGAGTACGCGCCACCACCGAATCCAACTCCGACCATCTAGTGGCTGCCAAT  
CCGTTTCGATGACAGCTACAACACCATTTCTATAAACCCTGCCTTCATCTAATCCATATCTTGGCCCTG  
GGTATCCTGGCTTTGGAGGCTACAGCACATTCAGAATGCCACCCACGTCCTCCAAGAATGTCTTCTCC  
CTACTGTGGTCTTACTCACTCAGGAATCAGCCCCACCCATTTCTCAGAATCCGTTGGGCATGGGCTTT  
AACCGGCTCATGCTTTAACTTTGGGCGCATGATAATTCGAATTTGGAAACCCACCTTATAATAATG  
TACTGACTCAGGACATTAACATGCCCGGTGACATTTTAGACAAGGCTCTGCTGAAAACCTCAGTCAGAT  
TCCCCCGCAGAATGTTGGCCAAGTGTCTAACCTGACCTCGCATCTAATTTGGCCCTGGGAATAATTCA  
AATTTTACCTCTCCGTTAGAAACGAATCATTGTTTATCCACCCCAACCGGTTTGGCCAAGCAAAG  
CTCCACTTCCAAACAAGACTTCACTCAAGGGCAACCAAAACCCGAATCAGAATTCGTCCTCACTACCC  
ACCTCACCTAAATATGGAGGATCCAGTCAATCAGAGTAACGTGAGTTAAAAATGTCAACAGAAACAAC  
GTTGTCCAAGAGAACAGCCGTTCCGGCAGCGCAGAGGCCACCAACCATGCGAATGGGACCCAGAACA  
AGCCCCGGCAGCCAGGGGCGCAGCTGACCTGTGCACCCCGACAAAAGCCGAAGTTCTCCTGCTCC  
CAGCCGGCATGGCCATTCTCCTCTGACCTGTGTACCCGTGCGGGATTTGTACAATGAAGTGAATGAC  
GATCAGGACGCCATTCTGTGTGAAGCCTTTGTGAGAAGTGGTTTCATCGCATCTGCACTGGAATGACCG  
AAACAGCTACGGGCTCCTGACAGCGGAAGCATCCGCAGTGTGGGGCTGTGACACGTGCATGGCTGACAA  
GGATGTCCAGCTCATGCGCACTAGAGAGGCCTTTGGTCCACCTGCCGTGGGCGGCGATGCC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR223833 representing NM\_028116  
 Red=Cloning site Green=Tags(s)


MSAEQDKPEIALKRVVGGDSGLDGLGGPNLQLGSPDKKKRKANTQGSFPPLSEYAPPPNPNSDHLVAAN  
 PFDDSYNTISYKPLPSSNPYLPGYPFGGGYSTFRMPPHVPPRMSSPYCGPYSLRNQPHFPQNPLGMGF  
 NRPFAFNFGPHDNSNFGNPPYNNVLTQDINMPGQHFQGSANF SQIPPQNVGQVSNPDLASNFAPGNNS  
 NFTSPLTNHSFIPPPNAFGQAKAPLPKQDFTQGATKTPNQSSSTHPPHLNMDPVPVQSNVELKNVNRNN  
 VVQENSRSGSAEATNNHANGTQNKPRQPRGAADLCTPDKSRKFSLLPSRHGHSSSDPVYPCGICTNEVND  
 DQDAILCEASCQKWFHRICTGMTETAYGLLLTAEASAVWGCDTCMADKDVQLMRTREAFGPPAVGGDA

SGPTRRPLEQKLI SEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9044\\_a02.zip](https://cdn.origene.com/chromatograms/mm9044_a02.zip)

**Restriction Sites:** SgfI-RsrII

**Cloning Scheme:**



**ACCN:** NM\_028116

**ORF Size:** 1251 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_028116.3](#)

**RefSeq Size:** 1745 bp

**RefSeq ORF:** 1254 bp

**Locus ID:** 72135


**UniProt ID:** [Q9D0P5](#)

**Cytogenetics:** 9 D

**MW:** 45.5 kDa

**Gene Summary:** Involved in signal transduction through the Wnt pathway.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR223833