

## Product datasheet for **MR223931L3V**

### Cstf2 (NM\_133196) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cstf2 (NM_133196) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cstf2
Synonyms:	C630034J23Rik; Cstf64; mir-3112
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_133196
ORF Size:	1740 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223931).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_133196.5</a> , <a href="#">NP_573459.1</a>
RefSeq Size:	4352 bp
RefSeq ORF:	1743 bp
Locus ID:	108062
UniProt ID:	<a href="#">Q8BIQ5</a>
Cytogenetics:	X E3
Gene Summary:	One of the multiple factors required for polyadenylation and 3'-end cleavage of mammalian pre-mRNAs. This subunit is directly involved in the binding to pre-mRNAs (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)