

Product datasheet for **MR223935L3V**

Nxn1 (NM_145598) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nxn1 (NM_145598) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nxn1
Synonyms:	A930031O08Rik; RdCVF; Txnl6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145598
ORF Size:	651 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223935).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145598.2 , NP_663573.1
RefSeq Size:	2785 bp
RefSeq ORF:	654 bp
Locus ID:	234404
UniProt ID:	Q8VC33
Cytogenetics:	8 B3.3
Gene Summary:	May play a role in cone cell viability, slowing down cone degeneration, does not seem to play a role in degenerating rods.[UniProtKB/Swiss-Prot Function]


[View online »](#)