

Product datasheet for **MR223941L3V**

Wisp3 (NM_001127376) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Wisp3 (NM_001127376) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Wisp3
Synonyms:	CCN6; ENSMUSG00000062074; Gm735; WISP-3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001127376
ORF Size:	1062 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223941).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001127376.1 , NP_001120848.1
RefSeq Size:	1418 bp
RefSeq ORF:	1065 bp
Locus ID:	327743
UniProt ID:	D3Z5L9
Cytogenetics:	10 B1
Gene Summary:	Plays a role in mitochondrial electron transport and mitochondrial respiration. [UniProtKB/Swiss-Prot Function]



[View online »](#)