

Product datasheet for **MR223966L3V**

Lpp (NM_001145952) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lpp (NM_001145952) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lpp
Synonyms:	9430020K16Rik; AA959454; AU024130; B130055L10Rik; C79715; D630048H16
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001145952
ORF Size:	1839 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223966).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001145952.1 , NP_001139424.1
RefSeq Size:	15462 bp
RefSeq ORF:	1842 bp
Locus ID:	210126
UniProt ID:	Q8BFW7
Cytogenetics:	16 B1



[View online »](#)

Gene Summary:

May play a structural role at sites of cell adhesion in maintaining cell shape and motility. In addition to these structural functions, it may also be implicated in signaling events and activation of gene transcription. May be involved in signal transduction from cell adhesion sites to the nucleus allowing successful integration of signals arising from soluble factors and cell-cell adhesion sites. Also suggested to serve as a scaffold protein upon which distinct protein complexes are assembled in the cytoplasm and in the nucleus (By similarity).
[UniProtKB/Swiss-Prot Function]