

## Product datasheet for **MR223993L3V**

### Jade1 (NM\_172303) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Jade1 (NM_172303) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Jade1
Synonyms:	AU041499; D530048A03Rik; mKIAA1807; Phf17
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_172303
ORF Size:	2502 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR223993).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_172303.4</a> , <a href="#">NP_758507.3</a>
RefSeq Size:	5584 bp
RefSeq ORF:	2505 bp
Locus ID:	269424
UniProt ID:	<a href="#">Q6ZPI0</a>
Cytogenetics:	3 B



[View online »](#)

**Gene Summary:**

Component of the HBO1 complex which has a histone H4-specific acetyltransferase activity, a reduced activity toward histone H3 and is responsible for the bulk of histone H4 acetylation in vivo. Transcriptional coactivator, it may also promote acetylation of nucleosomal histone H4 by KAT5. Promotes apoptosis. May act as a renal tumor suppressor. Negatively regulates canonical Wnt signaling; at least in part, cooperates with NPHP4 in this function (By similarity).[UniProtKB/Swiss-Prot Function]