

Product datasheet for MR224006L3

Usf1 (NM_009480) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Usf1 (NM_009480) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Usf1
Synonyms:	bHLHb1; bHLHb11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224006).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_009480
ORF Size:	930 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_009480.2 , NP_033506.1
RefSeq Size:	1822 bp
RefSeq ORF:	933 bp
Locus ID:	22278
UniProt ID:	Q61069
Cytogenetics:	1 79.4 cM
Gene Summary:	This protein encoded by this gene is a member of the basic-Helix-Hoop-Helix-Leucine zipper (bHLH-LZ) family and encodes a protein that can act as a transcription factor. Studies indicate that the basic region interacts with DNA at E-Box motifs, while the helix-loop-helix and leucine zipper domains are involved in dimerization with different partners. This protein is involved in a wide array of biological pathways, including cell cycle regulation, immune response, and responses to ultraviolet radiation. Mice lacking most of the coding exons of this gene often lacked both whiskers and nasal fur, and were prone to epileptic seizures, while mice lacking both this gene and another family member, Usf2, displayed embryonic lethality (PMID:9520440). Mutations in the human ortholog of this gene have been associated with Familial Combined Hyperlipidemia (FCHL) in humans. Pseudogenes of this gene are found on chromosome 11 and the X chromosome. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Mar 2015]

Product images:



Circular map for MR224006L3