

Product datasheet for **MR224007L3V**

Hnrnpul1 (NM_178089) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hnrnpul1 (NM_178089) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hnrnpul1
Synonyms:	E1B-AP5; E1BAP5; E130317O14Rik; Hnrnpul; Hnrpul1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_178089
ORF Size:	1020 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224007).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_178089.2
RefSeq Size:	1233 bp
RefSeq ORF:	1023 bp
Locus ID:	232989
UniProt ID:	Q8VDM6
Cytogenetics:	7 A3



[View online »](#)

Gene Summary:

Acts as a basic transcriptional regulator. Represses basic transcription driven by several virus and cellular promoters. When associated with BRD7, activates transcription of glucocorticoid-responsive promoter in the absence of ligand-stimulation. Plays also a role in mRNA processing and transport. Binds avidly to poly(G) and poly(C) RNA homopolymers in vitro (By similarity).[UniProtKB/Swiss-Prot Function]