

Product datasheet for **MR224012L3V**

Dydc1 (NM_027094) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dydc1 (NM_027094) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dydc1
Synonyms:	1700029M23Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_027094
ORF Size:	525 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224012).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_027094.1 , NP_081370.1
RefSeq Size:	995 bp
RefSeq ORF:	528 bp
Locus ID:	69496
UniProt ID:	Q9D9T0
Cytogenetics:	14
Gene Summary:	Plays a crucial role during acrosome biogenesis.[UniProtKB/Swiss-Prot Function]



[View online »](#)