

## Product datasheet for **MR224045L3V**

### Sema6d (NM\_199240) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sema6d (NM_199240) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sema6d
Synonyms:	1110067B02Rik; AA409156; D330011G23; mKIAA1479
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_199240
ORF Size:	3162 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224045).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_199240.2</a>
RefSeq Size:	6470 bp
RefSeq ORF:	3165 bp
Locus ID:	214968
UniProt ID:	<a href="#">Q76KF0</a>
Cytogenetics:	2 F1



[View online »](#)

**Gene Summary:**

Shows growth cone collapsing activity on dorsal root ganglion (DRG) neurons in vitro. May be a stop signal for the DRG neurons in their target areas, and possibly also for other neurons. May also be involved in the maintenance and remodeling of neuronal connections (By similarity).[UniProtKB/Swiss-Prot Function]