

Product datasheet for **MR224056**

Gorab (NM_178883) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gorab (NM_178883) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gorab
Synonyms:	NTKL-BP1; SCYL1-BP1; Scyl1bp1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR224056 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTCAGGATTGGGCGGGCTTCTCTGAGGAAGAGCTGAGGAGACTAAAGCAGAATAAAGATCCATTTG
AACCTCAGCGTCGAATTCCTGTGAAGAAACTCGACAACAGCTTCAGCGTGAAAAAGCCCTTCTAGAGCA
GAGCCAAAAGCTTGGCCTCCAAGATGGGTGAGCCTCATTGCTCCAGAGCAGCTGCTCTGCACCCAAA
CAGAGAGCTAACAGTCAAAGCCACGCTCTCCTCCCTGTGGCCCCAGTCCCCTCACACCCACCTCTT
CCTCTGGCGATGGAAGCTACCTGGTGTGGAAGTCAGCCCCAAGAACCAGGACTTGAGAATTCCCACCA
TGGTCACAAAAGTGCCGAGGTTCGAGCTCAAAGCCAGATTGCAAAGTGAGAAAAAGAAATGGAATTG
CAAGAGAAATCTCGTTGGGAAGTCTCCAACAAGAACAGCGACTAATGGAAGAGAAAAATAAACGTAAGA
AAGCCCTTTTGGCTCAAGCCATTGCGGAAAGATCTAAGAAGACCCAAGCAGAGACCATAAACTAAAGAG
GATCCAGAAGGAGTTGCAGGCCTTAGACGACATGGTGTGAGCTGACATTGGCATCCTCAGGAACAGAATC
GACCAGGCCAGCCTGGAGTATTCCTACGCACGGAAGCGCTTTGATAGAGCTGAAGCTGAGTACATTACAG
CAAAGCTGGATTTACAGCGCAAAACCGAGACGAAAGAGCAGCTCACCGAACACCTCTGCACCATCATCCA
GCAGAACGAGCTCCGAAAGGCCAAGAAGCTGGAGGAGCTGATGCAGCAGCTGGACGTGCAGGCTGATGAG
GAGGCTTTGCAACTCGAGGTGGAGTTGAACAGTTGCTTCGGGAACAGGAAGCAGAGGCAGCCAAACAAA
TGGCTTCTGTAGAGAGGCTGTGTCCGCCTGATGGCGAGAGTGTGTCTTCAGAACTGGCTGAAGAGAACAA
TGAGCCTCAAAGCAAGCTCCTTCCCCAGAGACAGACAAGCCAGGCAATGTTGCAGTAGCTCTCCCAT
AGACTAGACTGCCAGATCCAGGAGCCAAGAACTTTTCAGCTGCTGTGGCCACA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR224056 protein sequence
Red=Cloning site Green=Tags(s)



MAQDWAGFSEELRLKQNKDPFEPQRRIPVKKTRQQLQREKALLEQSQKLGLDQGSASLLPEQLLSAPK
QRANSQKPRSPSPVAPSPLTPTSSSGDGKLPVGVSQPQEPGLENSHHGKSAEVRAPKPDCKVEKKKME
QEKSRWEVLQQEQRLMEEKNKRRKALLAQIAERSKKTQAEITIKLKRIQKELQALDMVSADIGILRNRI
DQASLEYSYARKRFRDRAEAAYITAKL DLQRKTETKEQLTEHLCTIIQQNELRKAKKLEELMQQLDVQADE
EALQLEVEVEQLLRQEAEAAKQMASVERLCPDGEVSSELAENNEPQKQAPSPETDKPGKCCSSSPH
RLDCPDGPAGKANSAEAT

TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM 178883

ORF Size: 1104 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_178883.6, NP_849214.2

RefSeq Size: 2509 bp

RefSeq ORF: 1107 bp


Locus ID: 98376

UniProt ID: Q8BRM2

Cytogenetics: 1 H2.1

MW: 41.5 kDa

Product images:



Circular map for MR224056