

Product datasheet for **MR224065L3V**

Herc6 (NM_025992) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Herc6 (NM_025992) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Herc6
Synonyms:	1700121D12Rik; 2510038N07Rik; 4930427L17Rik; AI451296; CEB1; Herc5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025992
ORF Size:	3009 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224065).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025992.2 , NP_080268.1
RefSeq Size:	5279 bp
RefSeq ORF:	3012 bp
Locus ID:	67138
UniProt ID:	F2Z461
Cytogenetics:	6
Gene Summary:	Major E3 ligase for ISG15 conjugation. Acts as a positive regulator of innate antiviral response in cells induced by interferon.[UniProtKB/Swiss-Prot Function]



[View online »](#)