

Product datasheet for **MR224071L3V**

Mc3r (NM_008561) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mc3r (NM_008561) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mc3r
Synonyms:	MC3-R
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008561
ORF Size:	972 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224071).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008561.3 , NP_032587.1
RefSeq Size:	2623 bp
RefSeq ORF:	972 bp
Locus ID:	17201
UniProt ID:	P33033
Cytogenetics:	2 94.8 cM



[View online »](#)

Gene Summary:

This gene encodes a member of the melanocortin receptor family. Melanocortin receptors are transmembrane G-protein coupled receptors, which respond to small peptide hormones and exhibit diverse functions and tissue type localization. As part of the central nervous melanocortin system, the encoded protein is competitively bound by either melanocyte stimulating hormone or agouti-related protein to regulate energy homeostasis and adaptation to food restriction. Disruption of this gene results in an increased ratio of weight gain to food intake, increased fat mass, and decreased lean mass, without having a large effect on insulin sensitivity or glucose metabolism. [provided by RefSeq, Dec 2012]