

## Product datasheet for **MR224083L3V**

### **Hmbox1 (NM\_177338) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Hmbox1 (NM_177338) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hmbox1
Synonyms:	AI451877; AI604847; F830020C16Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_177338
ORF Size:	1257 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224083).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_177338.5</a>
RefSeq Size:	2931 bp
RefSeq ORF:	1260 bp
Locus ID:	219150
UniProt ID:	<a href="#">Q8BJA3</a>
Cytogenetics:	14 D1



[View online »](#)

**Gene Summary:**

Binds directly to 5'-TTAGGG-3' repeats in telomeric DNA (By similarity). Associates with the telomerase complex at sites of active telomere processing and positively regulates telomere elongation (By similarity). Important for TERT binding to chromatin, indicating a role in recruitment of the telomerase complex to telomeres (PubMed:23685356). Also plays a role in the alternative lengthening of telomeres (ALT) pathway in telomerase-negative cells where it promotes formation and/or maintenance of ALT-associated promyelocytic leukemia bodies (APBs) (By similarity). Enhances formation of telomere C-circles in ALT cells, suggesting a possible role in telomere recombination (By similarity). Might also be involved in the DNA damage response at telomeres (By similarity).[UniProtKB/Swiss-Prot Function]