

## Product datasheet for **MR224159L3V**

### Tcf12 (NM\_011544) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tcf12 (NM_011544) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tcf12
Synonyms:	A130037E08Rik; ALF1; bHLHb20; HEB; HEBAIt; HTF-4; HTF4; ME1; REB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011544
ORF Size:	2121 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224159).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011544.3</a> , <a href="#">NP_035674.2</a>
RefSeq Size:	4707 bp
RefSeq ORF:	2121 bp
Locus ID:	21406
UniProt ID:	<a href="#">Q61286</a>
Cytogenetics:	9 39.85 cM
Gene Summary:	Transcriptional regulator. Involved in the initiation of neuronal differentiation. Activates transcription by binding to the E box (5'-CANNTG-3').[UniProtKB/Swiss-Prot Function]



[View online »](#)