

Product datasheet for **MR224161L3V**

Slc20a2 (NM_011394) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc20a2 (NM_011394) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc20a2
Synonyms:	MolPit2; Pit-2; Pit2; Ram-1; Ram1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011394
ORF Size:	1968 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224161).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011394.3 , NP_035524.2
RefSeq Size:	3658 bp
RefSeq ORF:	1971 bp
Locus ID:	20516
UniProt ID:	Q80UP8
Cytogenetics:	8 11.42 cM



[View online »](#)

Gene Summary:

Sodium-phosphate symporter which seems to play a fundamental housekeeping role in phosphate transport by absorbing phosphate from interstitial fluid for normal cellular functions such as cellular metabolism, signal transduction, and nucleic acid and lipid synthesis. In vitro, sodium-dependent phosphate uptake is not significantly affected by acidic and alkaline conditions, however sodium-independent phosphate uptake occurs at acidic conditions. May play a role in extracellular matrix, cartilage and vascular calcification. Functions as a retroviral receptor (By similarity).[UniProtKB/Swiss-Prot Function]