

Product datasheet for **MR224167L4V**

Adarb1 (NM_001024837) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adarb1 (NM_001024837) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Adarb1
Synonyms:	1700057H01Rik; AD; Adar2; AW124433; AW558573; BB220382; D10Bwg0447e; RED; Red1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001024837
ORF Size:	2133 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224167).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001024837.2 , NP_001020008.1
RefSeq Size:	6602 bp
RefSeq ORF:	2136 bp
Locus ID:	110532
UniProt ID:	Q91ZS8
Cytogenetics:	10 39.72 cM



[View online »](#)

Gene Summary:

This gene encodes a double-stranded-RNA-specific adenosine deaminase that is involved in editing pre-mRNAs by site-specific conversion of adenosine (A) to inosine (I). Substrates for this enzyme include ionotropic glutamate receptors (GluR2-6) and serotonin receptor (5HT2C). Studies in rodents have shown that this protein can modify its own pre-mRNA by A->I editing to create a novel acceptor splice site, alternative splicing to which results in down regulation of its protein expression. Additional splicing events result in transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]