

## Product datasheet for **MR224174L4V**

### **Ilk (NM\_001161724) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ilk (NM_001161724) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ilk
Synonyms:	AA511515; ESTM24
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001161724
ORF Size:	1356 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224174).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001161724.1</a> , <a href="#">NP_001155196.1</a>
RefSeq Size:	1797 bp
RefSeq ORF:	1359 bp
Locus ID:	16202
UniProt ID:	<a href="#">O55222</a>
Cytogenetics:	7 E3



[View online »](#)

**Gene Summary:**

Receptor-proximal protein kinase regulating integrin-mediated signal transduction. May act as a mediator of inside-out integrin signaling. Focal adhesion protein part of the complex ILK-PINCH. This complex is considered to be one of the convergence points of integrin- and growth factor-signaling pathway. Could be implicated in mediating cell architecture, adhesion to integrin substrates and anchorage-dependent growth in epithelial cells. Phosphorylates beta-1 and beta-3 integrin subunit on serine and threonine residues, but also AKT1 and GSK3B.[UniProtKB/Swiss-Prot Function]