

Product datasheet for **MR224176L4V**

A1cf (NM_001081074) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	A1cf (NM_001081074) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	A1cf
Synonyms:	1810073H04Rik; Acf; ACF64; ACF65; ASP; MCM; mer; MerCreMer; Tg(Myh6-cre/Esr1*)1Jmk
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001081074
ORF Size:	1785 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224176).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001081074.1
RefSeq Size:	3828 bp
RefSeq ORF:	1788 bp
Locus ID:	69865
UniProt ID:	Q5YD48
Cytogenetics:	19 C1



[View online »](#)

Gene Summary:

Essential component of the apolipoprotein B mRNA editing enzyme complex which is responsible for the postranscriptional editing of a CAA codon for Gln to a UAA codon for stop in APOB mRNA. Binds to APOB mRNA and is probably responsible for docking the catalytic subunit, APOBEC1, to the mRNA to allow it to deaminate its target cytosine. The complex also seems to protect the edited APOB mRNA from nonsense-mediated decay (By similarity).
[UniProtKB/Swiss-Prot Function]