

Product datasheet for **MR224182L4V**

Tcf7l2 (NM_001142922) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Tcf7l2 (NM_001142922) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Tcf7l2 |
| Synonyms: | Tcf-4; Tcf4; TCF4B; TCF4E |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001142922 |
| ORF Size: | 1794 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR224182). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001142922.1 |
| RefSeq Size: | 4138 bp |
| RefSeq ORF: | 1797 bp |
| Locus ID: | 21416 |
| Cytogenetics: | 19 51.59 cM |



[View online »](#)

Gene Summary:

Participates in the Wnt signaling pathway and modulates MYC expression by binding to its promoter in a sequence-specific manner. Acts as repressor in the absence of CTNNB1, and as activator in its presence. Activates transcription from promoters with several copies of the Tcf motif CCTTTGATC in the presence of CTNNB1. TLE1, TLE2, TLE3 and TLE4 repress transactivation mediated by TCF7L2/TCF4 and CTNNB1. Expression of dominant-negative mutants results in cell-cycle arrest in G1 (By similarity). Necessary for the maintenance of the epithelial stem-cell compartment of the small intestine.[UniProtKB/Swiss-Prot Function]