

## Product datasheet for MR224202

### Mblac2 (NM\_028372) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Mblac2 (NM_028372) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Mblac2
Synonyms:	2900024O10Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR224202 representing NM_028372 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCGGCGCTCGAGTGGTATGCCACAAGTCTCTGGCGATGGCATCTTCTGGATTCAAGAACGGTTCT  
ACGAGTCAGGCAACCGGCCAACATTTGGCTGGTGC GCGGCTCGGAGCAGGACGTGGTGATAGACTGG  
CCTGGGGCTGCGCAGCCTTCCGGAGTACCTCTACTCCTCTGGCCTTTCGAGGATTGCGGGTCTAAAGAG  
GATGCAGGACGCCGCCACTGCTGGCCGTGGCCACCCACGTGCACTTCGATCACTCCGGCGGTCTCTATC  
AGTTCGACCAGGTGGCTGTGCACCGCGCCGAGGCCGAGGCCCTGGCTCGCGGGGACAACCTTGAGACCGT  
GACCTGGCTTTCTGACAGCGAGGTGGTGC GGGCGCCAGCCCTGGCTGGAGGGCCAGGCAGTTCCGAGTG  
CAGGCAGTGCAACCCACCCTCATTCTGCAGGATGGGGATGTGATCAACCTTGCGCAGACAGACTCACTG  
TTATGCACATGCCAGGTCACTCCAGAGGCAGTATTTGCTTACATGATAAAGACCGTAAGGTGCTCTTCAG  
TGGTGACGTGCTGATGATGGATCACTGATTGACTGGCTCCCATACAGCAGGATCAGTGACTATGTTGGA  
ACTTGTGAGCGTCTCATAGAGTTGGTGGACCGAGGCCCTTGTGGAGAAGGTGCTTCTGGGCACTTCAATA  
CCTTTGGTGCCGAAAGGCTGTTTCGCTTAGCTTAACTATATTTCAAAGCTGGGATATGCCACAAAGT  
GTCTACTTTTGCCATGAGGTCTCTTGAAGTTTAGCTCTCCGTGAACAAATCTAGGACCTCACCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR224202 representing NM\_028372  
 Red=Cloning site Green=Tags(s)

MSALEWYAHKSLGDGIFWIQERFYESGNRANIWLVRGSEQDVVIDTGLGLRSLPEYLYSSGLLQDCGSKE  
 DAGRRPLLAVATHVHFDHSGGLYQFDQVAVHRAEAEALARGDNFETVTWLS DSEVVVRAPSPGWRARQFRV  
 QAVQPTLILQGDVINLGDRLTVMHMPGHSRGSICLHDKDRKVLFSGDVVYDGSLLDWLPYSRISDYVG  
 TCERLIELVDRGLVEKVLPGHFNTFGAERLFRSLASNYISKAGICHKVSTFAMRSLASLALRVTNPRTSP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9079\\_a08.zip](https://cdn.origene.com/chromatograms/mm9079_a08.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_028372

**ORF Size:** 837 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_028372.1](#), [NP\\_082648.1](#)

**RefSeq Size:** 3881 bp

**RefSeq ORF:** 840 bp


**Locus ID:** 72852

**UniProt ID:** [Q8BL86](#)

**Cytogenetics:** 13 C3

**MW:** 31.2 kDa

**Product images:**



Circular map for MR224202