

## Product datasheet for **MR224221L3V**

### Ilf3 (NM\_010561) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ilf3 (NM_010561) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ilf3
Synonyms:	MBII-26; MPHOSPH4; NF9; NF90; NFAR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010561
ORF Size:	2733 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224221).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_010561.2</a> , <a href="#">NP_034691.2</a>
RefSeq Size:	3457 bp
RefSeq ORF:	2736 bp
Locus ID:	16201
UniProt ID:	<a href="#">Q9Z1X4</a>
Cytogenetics:	9 7.78 cM



[View online »](#)

**Gene Summary:**

The protein encoded by this gene contains two double-stranded RNA binding domains and functions in the post-transcriptional regulation of gene expression. It is a component of an RNA-protein complex that may be involved in mediating the export of messenger RNAs. Alternative splicing results in multiple transcript variants encoding distinct isoforms. These isoforms are grouped into two categories, NFAR-1 or NFAR-2, based on variation at the C-terminus. [provided by RefSeq, Mar 2013]