

Product datasheet for **MR224229**

Ilf3 (NM_001042709) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ilf3 (NM_001042709) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ilf3
Synonyms:	MBII-26; MPHOSPH4; NF9; NF90; NFAR
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR224229 representing NM_001042709
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCGTCCCATGAGAATTTTTGTGAATGATGATCGCCACGTGATGGCAAAGCATTCTTCAGTGTACCCAA
 CACAAGAGGAGCTGGAGGCTGTACAGAACATGGTGTCCCATACTGAGCGGGCCCTGAAGGCTGTCTCTGA
 CTGGATTGATGAGCAGGAGAAAAGCAACAGCGAGCTCTCTGAGGCAGAAAATATGGACACACCCCCAGAC
 GATGAGAGCAAAGAAGGGGCTGGGGAACAGAAGGCGGAACACATGACTAGGACCCCTGAGGGGCGTGATGC
 GGGTCGGCCTGGTAGCCAAGGGTCTTCTGCTCAAGGGGGACTTGGATCTGGAGCTGGTTCTGCTCTGTAA
 GGAGAAGCCACAACCGCCCTTCTGGACAAGGTGGCTGACAACCTGGCCATCCAGCTCACTACTGTAAACA
 GAAGACAAGTATGAAATACTCCAGTCTGTGGATGATGCTGCGATTGTGATAAAAAACACAAAAGAGCCCC
 CCTTGTCTTGGACATCCATCTGACCTCCCTGTTGTGACAGAGAAGAAATGGAGAAAGTATTAGCTGGAGA
 AACGCTATCAGTCAACGATCCCCGGACGTTCTGGACAGGCAGAAATGCCTTGTGCCTTGGCGTCCCTC
 CGACACGCCAAGTGGTTCAGGCCAGAGCCAATGGACTGAAGTCATGTGTATTGTCATCCGTGTCTTAA
 GGGACTTGTGTACCCGAGTGCCACCTGGGGTCCCTCAGAGGATGGCCTCTGGAGCTGTGTGTGAGAA
 GTCCATCGGCACTGCCAATAGGCCAATGGGTGCTGGTGAAGCCCTGCGGAGAGTGTGGAGTGCCTGGCA
 TCCGGCATCGTAATGCCAGATGGTCTGGCATTATGACCCTTGTGAAAAAGAAGCCACTGATGCTATTG
 GGCATCTAGACAGACAGCAACGGGAAGATACACAGAGTGCAGCAGCATGCTCTGCGGCTTGTGCCTT
 TGGTCAACTCCATAAAGTACTGGGAATGGACCCCTGCCTTCCAAAATGCCCAAGAAACCAAGAACGAG
 AACCCGGTGGACTACACTGTTCAAATTCCTCCAGCACCACCTATGCTATCACACCATGAAACGCCCTA
 TGAAGAGGATGGGGAGGAGAAGTCTCCAGCAAGAAGAAAAAGAAAGATCCAGAAGAAAGAGGAGAAGGC
 TGATCCTCCTCAAGCTATGAATGCCCTGATGAGGTTAAATCAGCTGAAGCCAGGGCTGCAGTACAAGCTG
 ATCTCCAGACAGGCCCTGTTTATGCTCCCATCTTACCATGTCTGTGGAGGTAGACGGCAGTAACTTCG
 AGGCCTCGGGGCCATCTAAAAAGACTGCCAAGCTTATGTAGCTGTGAAGGTGTTACAGGACATGGGCTT
 GCCAACAGGCGCTGAAGGCAGAGACTCCAGCAAGGGGGAAGACTCCGCTGAGGAGTCAGATGGGAAGCCA
 GCAATAGTGGCCCCACCCCTGTGGTGAAGCTGTCTCCAACCCAGTCTGTCTTCCCTTCAGATGCCA
 CTACTGAGCAGGGACCGATTTTACTAAGCATGGCAAGAACCCTGTTATGGAGCTTAATGAGAAGAGACG
 TGGCCTCAAATATGAGCTCATTTCTGAGACGGGGGCAGCCACGACAAAAGGTTTGTATGGAGGTTGAG
 GTGGACGGACAGAAGTTTCAAGGTGCTGGTTCAAACAAAAGGTGGCAAAGGCTTATGCTGCACTTGCCG
 CATTAGAAAAACTTTTCCCTGATACCCCTTGTCTTGAAGCCAACAAAAGAAAAGGACCCCACTACC
 TGTCCGAGGTGGACCAATTTGCTGCCAAGCCACACAACCCTGGTTTTTGGCATGGGAGGCCCATGCAT
 AATGAAGTGCCGCCACCTCCTAACATCCGAGGTCGGGGCCGAGGAGGTAACATCCGAGGGCGAGGACGGG
 GGCGAGGATTTGGTGGCGCAACCATGGAGGAGGCTACATGAATGCTGGTGTGGATATGGAAGCTATGG
 GTACAGCAGCAATTCGGCCACAGCAGGCTACAGTGACTTTTTACAGACTGCTACGGCTATCATGATTTT
 GGGGCTTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR224229 representing NM_001042709
Red=Cloning site Green=Tags(s)


MRP M R I F V N D D R H V M A K H S S V Y P T Q E E L E A V Q N M V S H T E R A L K A V S D W I D E Q E K G N S E L S E A E N M D T P P D
 D E S K E G A G E Q K A E H M T R T L R G V M R V G L V A K G L L L K G D L D E L V L L C K E K P T T A L L D K V A D N L A I Q L T T V T
 E D K Y E I L Q S V D D A A I V I K N T K E P P L S L T I H L T S P V V R E E M E K V L A G E T L S V N D P P D V L D R Q K C L A A L A S L
 R H A K W F Q A R A N G L K S C V I V I R V L R D L C T R V P T W G P L R G W P L E L L C E K S I G T A N R P M G A G E A L R R V L E C L A
 S G I V M P D G S G I Y D P C E K E A T D A I G H L D R Q Q R E D I T Q S A Q H A L R L A A F G Q L H K V L G M D P L P S K M P K P K N E
 N P V D Y T V Q I P P S T T Y A I T P M K R P M E E D G E E K S P S K K K K I Q K K E E K A D P P Q A M N A L M R L N Q L K P G L Q Y K L
 I S Q T G P V H A P I F T M S V E V D G S N F E A S G S K K T A K L H V A V K V L Q D M G L P T G A E G R D S S K G E D S A E E S D G K P
 A I V A P P P V V E A V S N P S S V F P S D A T T E Q G P I L T K H G K N P V M E L N E K R R G L K Y E L I S E T G G S H D K R F V M E V E
 V D G Q K F Q G A G S N K K V A K A Y A A L A E K L F P D T P L A L E A N K K K R T P V P V R G G P K F A A K P H N P G F G M G G P M H
 N E V P P P N I R G R G R G N I R G R G R G F G G A N H G G G Y M N A G A G Y G S Y G Y S S N S A T A G Y S D F F T D C Y G Y H D F
 G A S

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001042709

ORF Size: 2109 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001042709.2](#), [NP_001036174.1](#)

RefSeq Size: 3606 bp

RefSeq ORF: 2112 bp


Locus ID: 16201

Cytogenetics: 9 7.78 cM

MW: 76.8 kDa

Gene Summary: The protein encoded by this gene contains two double-stranded RNA binding domains and functions in the post-transcriptional regulation of gene expression. It is a component of an RNA-protein complex that may be involved in mediating the export of messenger RNAs. Alternative splicing results in multiple transcript variants encoding distinct isoforms. These isoforms are grouped into two categories, NFAR-1 or NFAR-2, based on variation at the C-terminus. [provided by RefSeq, Mar 2013]

Product images:



Circular map for MR224229