

Product datasheet for **MR224232L4V**

Dok7 (NM_172708) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dok7 (NM_172708) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dok7
Synonyms:	A930013K19Rik; AW049091; Dok-7; EF-12; Oit5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172708
ORF Size:	1710 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224232).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172708.3
RefSeq Size:	2483 bp
RefSeq ORF:	1713 bp
Locus ID:	231134
Cytogenetics:	5 B2
Gene Summary:	Probable muscle-intrinsic activator of MUSK that plays an essential role in neuromuscular synaptogenesis. Acts in aneural activation of MUSK and subsequent acetylcholine receptor (AChR) clustering in myotubes. Induces autophosphorylation of MUSK.[UniProtKB/Swiss-Prot Function]



[View online »](#)