

## Product datasheet for **MR224234L3V**

### Gm7325 (Mymx) (NM\_001177470) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Gm7325 (Mymx) (NM_001177470) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mymx
Synonyms:	EG653016; Esgp; Gm7325; minion; myomerg
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001177470
ORF Size:	324 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224234).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001177470.1</a> , <a href="#">NP_001170941.1</a>
RefSeq Size:	852 bp
RefSeq ORF:	327 bp
Locus ID:	653016
UniProt ID:	<a href="#">Q2Q5T5</a>
Cytogenetics:	17 B3



[View online »](#)

**Gene Summary:**

Myoblast-specific protein that mediates myoblast fusion, an essential step for the formation of multi-nucleated muscle fibers (PubMed:28386024, PubMed:28569745, PubMed:28569755, PubMed:30197239). Involved in membrane fusion downstream of the lipid mixing step mediated by MYMK (PubMed:30197239). Acts by generating membrane stresses via its extracellular C-terminus, leading to drive fusion pore formation (PubMed:30197239). Acts independently of MYMK (PubMed:30197239). Involved in skeletal muscle regeneration in response to injury by mediating the fusion of satellite cells, a population of muscle stem cells, with injured myofibers (PubMed:29581287).[UniProtKB/Swiss-Prot Function]