

Product datasheet for **MR224272**

Klhdc1 (NM_178253) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Klhdc1 (NM_178253) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Klhdc1
Synonyms:	MGC141300; MGC141301
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR224272 representing NM_178253 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGCGGACCCGAGGCGTTCTGCGTGGCCGAGGAGCGCAGCGGCCACTGCGCAGTGGTGGACGGACACT
TCCTCTACGTGTGGGGGGCTACGTGTCTATTGAAGACAATGAGGTATTTGCCAAATGATGAAATGTG
GACCTACGATATTGATAGTGGGCTGTGAAAATGCACCTAATGGAAGGGAGCTCCCCCTCCATGTG
GGAAGCTGCGGGCTTGATCCACGGAAGGCTCTATGTCTTTGGAGGCTATGACGACAAAGGATACAGCA
ACCGGCTTTACTTTGTCAATTTACGAACAAGAGATGGAACCTACACCTGGGAAAAATCACCAAGTTTGA
TGGACAACCACCCACACCACGTGATAAACTTTCTGCTGGGTATATAAAGACAGATTAATATATTTTCGGT
GGTTATGGGTACAGGAGGCACAGTGAGCTTCAAGAATGTTTTGATGTTTCATGATGCATCATGGGAAGAAC
AGATATTCTGGGGGTGGCACAATGACGTCCATGTGTTTGACACAAAGACGAGGACGTGGTCTCAACCAGA
AATTAAGGTGGAGTTCACCACAGCCACGAGCCGCGCACTCATGTGCAGTTCTGGGAAATAAGGGTTAC
GTCTTTGGCGGACGTGTTTTGCAAACAGGATGAATGACCTACACTATCTAAACCTAGACACCTGGGTCT
GGTCTGGAAGGATTTCTGTTAATGGAGAAAGCCCCAACACCGGTCATGGCACACTTTAACAGCAATAAC
TGATGATAAGCTTTCTATTTGGTGGACTAAATGCAGATAATATCCACTAAGTGATGGCTGGATTTCAT
AATATAACAACAACTGTTGAAACAACCTAGACATTTACCGTACACCAGGCCAGGTTATGGCACACAG
CCTGTCTGGGAAAAGAAAATGAAATAATGGTATTGGTGGGAGCAAAGACAACCTTACTTTCTTGATAC
AGGTCACTGTAATGATTTATTGATCTTTCAAACACAGCCCTATTCACTGCTCAGGTCATGTCTTGACTGC
ATTGGGAAAAATGCTATCATTTTAAAAAGTCAAATATCTTTATTACCTCCTAACTTCTGCAGCAAGTCC
TCAAAAAAATAACATTTTGGACTGCAGCTAATTACCGAAAGGAGCAAAGAATCCGAAAGGAAGAGACAGA
AAATAATCAGCCGAGAGTCAGTAGCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR224272 representing NM_178253
 Red=Cloning site Green=Tags(s)


MADPQAFCAEERSGHCAVVDGHFLYVWGGYVSIEDNEVYLPNDEMWTYDIDSGLWKMHLMEGELPPSMS
 GSCGACIHGRLYVFGGYDDKGYSNRLYFVNLRTDGTYTWEKITKFDGQPPTPRDKLSCWVYKDRLIYFG
 GYGYRRHSELQECFDVHDASWEEQIFWGWHDVHVFDTKTRTWSQPEIKGGVPPQPRAAHSCAVLGNKGY
 VFGGRVLQTRMNDLHYLNLDTWVWSGRISVNGESPKHRSWHTLTAITDDKLFLFGGLNADNIPLSDGWIH
 NITTCNKWQLRHLPTYTRPRLWHTACLKENEIMVFGGSKDNLLFLDTGHCNDLLIFQTQPYSLLRCLDC
 IGKNAIILKSQISLLPPKLLQVLKKITFWTAANYRKEQIRKEETENNQPRVSSC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9070_g03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_178253

ORF Size: 1218 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_178253.5, NP_839984.1

RefSeq Size: 1949 bp

RefSeq ORF: 1221 bp


Locus ID: 271005

UniProt ID: Q80YG3

Cytogenetics: 12 C2

MW: 46.8 kDa

Product images:



Circular map for MR224272