

## Product datasheet for **MR224295L3V**

### Hif3a (NM\_016868) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hif3a (NM_016868) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hif3a
Synonyms:	bHLHe17; lpas; MOP7; NEPAS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016868
ORF Size:	1989 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224295).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_016868.3</a> , <a href="#">NP_058564.2</a>
RefSeq Size:	6947 bp
RefSeq ORF:	1989 bp
Locus ID:	53417
UniProt ID:	<a href="#">Q0VBL6</a>
Cytogenetics:	7 A2



[View online »](#)

**Gene Summary:**

Acts as a transcriptional regulator in adaptive response to low oxygen tension. Acts as a regulator of hypoxia-inducible gene expression (PubMed:9840812, PubMed:11734856, PubMed:21546903). Plays a role in the development of the cardiorespiratory system (PubMed:18070924).[UniProtKB/Swiss-Prot Function]