

Product datasheet for **MR224297L3V**

Klf13 (NM_021366) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Klf13 (NM_021366) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Klf13
Synonyms:	0610043C13Rik; 9430029L20Rik; Bteb3; FKLF-2; FKLF2; NSLP1; RFLAT-1; RFLAT1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_021366
ORF Size:	867 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224297).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021366.3 , NP_067341.2
RefSeq Size:	6396 bp
RefSeq ORF:	870 bp
Locus ID:	50794
UniProt ID:	Q9JIZ6
Cytogenetics:	7 C
Gene Summary:	Transcription factor that activates expression from GC-rich minimal promoter regions, including genes in the cells of the erythroid lineage.[UniProtKB/Swiss-Prot Function]



[View online »](#)