

Product datasheet for **MR224299L4V**

Eda (NM_001177938) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Eda (NM_001177938) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Eda |
| Synonyms: | Ed1; Eda-A1; Eda-A2; EDA1; HED; Ta; tabby; Tnlg7c; XLHED |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001177938 |
| ORF Size: | 1164 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR224299). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001177938.1 , NP_001171409.1 |
| RefSeq Size: | 4946 bp |
| RefSeq ORF: | 1167 bp |
| Locus ID: | 13607 |
| Cytogenetics: | X 43.59 cM |



[View online »](#)

Gene Summary:

Cytokine which is involved in epithelial-mesenchymal signaling during morphogenesis of ectodermal organs. Functions as a ligand activating the DEATH-domain containing receptors EDAR and EDA2R. Isoform TAA binds only to the receptor EDAR, while isoform TA-A2 binds exclusively to the receptor EDA2R (By similarity). May also play a role in cell adhesion (PubMed:10534613).[UniProtKB/Swiss-Prot Function]