

Product datasheet for **MR224307L4V**

Sap30 (NM_021788) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sap30 (NM_021788) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sap30
Synonyms:	30kDa
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021788
ORF Size:	663 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224307).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021788.2 , NP_068560.1
RefSeq Size:	1183 bp
RefSeq ORF:	663 bp
Locus ID:	60406
UniProt ID:	O88574
Cytogenetics:	8 29.85 cM



[View online »](#)

Gene Summary:

Involved in the functional recruitment of the Sin3-histone deacetylase complex (HDAC) to a specific subset of N-CoR corepressor complexes. Capable of transcription repression by N-CoR. Active in deacetylating core histone octamers (when in a complex) but inactive in deacetylating nucleosomal histones.[UniProtKB/Swiss-Prot Function]