

Product datasheet for **MR224334A1V**

Mouse Vmn1r200 (NM_134212) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Vmn1r200 (NM_134212) AAV Particle
Tag:	Myc-DDK
Symbol:	Vmn1r200
Synonyms:	V1rh3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR224334 representing NM_134212 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGATTTTAAAGTACATTAAAGAAATTATTTCCCTCCATGACAGTGCTTGGCACTCTGGGAAACATGT
TTGTTTTGTGAATTATATGTTCAAGTTGGTGGGAAGCCCTGAGAAGAAACCATACACCTCATTCTCAT
CCACTTGGCTTTTACAAACATCATAATCCTTCTTACAAAAGGATTGAAAAGACAATAGCAGCTTTTGGT
TTGAGAACTTCTGGATGACATAGGCTGTAAGAGCATTGTCTACCTGGAAAGGTATCCCGTGGCCTCT
CCATCTGTACCAGCAGTCTCCTCACTGTGGTCCAGGCCATCACCATCAGTCCAGAGCATCTGGATGGAG
GAGGCTCAGGACAAAGTCTGCATGGCACATCCTTCCATTCTTGCTATTCTTTTGGATACTCAATGGTTTA
ATAAGTATGAACCTAATCCATTCCATCACAGGTACAGGTTTTAATATATCACAGCTTGAGAATAGAGACA
GCTATTGCTATTTTATGCTAGAAAGTCAGGAAATAAATGGATTGTTCTCCCTCTCATGGTCATGAGAGA
TGCAGTGTTTCAGAGTGCCCTGGGAGGGGCCAGTGGCTACATGATATTCTTCTCCACAAGCACCACCAG
CATGTCCTCTACCTTCAGAAATCCAAGTTCTCTACAGAACTCCCCCTGAGCTGAGAGCTGCCAGAGTG
TCCTCCTCTGATGCTCTGTTTTGCTTTCTTCTACTGGACAGATTGTGCGTTTTCTCAATTTTAAAGTCT
CTCTTCAAGGGACAAAACCTTGATGGTAAATACTCAAGAATTTCTGGCCCTTGGTTATGCAACTTTTAGC
CCTTTTGTGCTGATTCACAGGGATGGGCTTCTGGTTGATTGGTGGCATGCTCAGATGGAGAAATTGAGAA
AATGTTACTCTGTTTATATGTCAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR224334 representing NM_134212 Red=Cloning site Green=Tags(s) MILKYIKEIIFPSMTVLGTLGNMFVFNMFSSWWSPEKKPIHLILIHLAFTNIIILLTKGLQKTIAAFGLRNFLDDIGCKSIYVLERSRGLSICTSLLTVVQAITISPRASGWRRRLRTKSAWHILPFLFFWILNGLISMNLIHSITGTGFNISQLENRDSYCYFMLESQEIKWIVLPLMVMRDAVFQSALGGASGYMIFLLHKHHQHVLYLQNSKLLYRTPPELRAAQSVLLMLCFAFFYWTDCAFSLILSLSSRDKTLMVNTQEFALALGYATFSPFVLIHRDGLLDVWWHAQMEKLRKCYSRLYVQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_134212
ORF Size:	936 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_134212.1, NP_598973.1</u>
RefSeq Size:	939 bp
RefSeq ORF:	939 bp
Locus ID:	171246
Cytogenetics:	13 A3.1
MW:	36.1 kDa