

## Product datasheet for **MR224352L4V**

### Gpr26 (NM\_173410) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Gpr26 (NM_173410) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gpr26
Synonyms:	9630036A11Rik; A730021D22
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_173410
ORF Size:	1011 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224352).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_173410.3</a> , <a href="#">NP_775586.1</a>
RefSeq Size:	2799 bp
RefSeq ORF:	1014 bp
Locus ID:	233919
UniProt ID:	<a href="#">Q8BZA7</a>
Cytogenetics:	7 F3
Gene Summary:	Orphan receptor. Displays a significant level of constitutive activity. Its effect is mediated by G(s)-alpha protein that stimulate adenylate cyclase, resulting in an elevation of intracellular cAMP (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)