

## Product datasheet for **MR224354**

### Arhgap27 (NM\_133715) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Arhgap27 (NM_133715) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Arhgap27
Synonyms:	2310069I04Rik; 5730442P18Rik; Gm21552; Sh3d20
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MR224354 representing NM\_133715  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGTGGATATGATTTCAAACCTGGTCAGAAGGCAGAGCCGGGCCCTGCGGGCACAGGTGGACGACCCCC  
 CGGAACCAAGTGTATGCCAATGTAGAGAGGCAACCTCGGGCCACGTCGCCTCGCTCAGCTGCAGCTCCTCC  
 TCGGCTCAGCCCAGTGTGGGAGACACACAGACACGGGCACTGGGCGCCCTACTACTACAACCCCGAC  
 ACGGGCGTGACCACCTGGGAGTCAACCTTTGAAACACCTGAAGGCACCACCAGCCCGCCACCTCCCGGG  
 CCTCTGTGGTAGCGGGGAGAGCCTGGAGACAGAGTGGGGCCAGTACTGGGATGAGGAGAGCCGAGGGT  
 CTTCTTACAACCCCTTGACCGGTGAGACGGCCTGGGAGGACGAGACCGAGGAGCTGGAGGAAGATCAC  
 CAGGAACAACCTGAAATGCAACCCAGCCTGAGCCCGCAAGCCCCGGGACGAGAGGCCCAACTCCTG  
 AAACAGACTACCCGAATTACTGGCCAGTACCCTGAAGAGGACTACTCTCTGTGGGCTCTTCAGTGA  
 TCCGGGCCCGCTTCTCCTTTGGTTGCACCCCTGGCTGGTCTGTCAAATTACCCAGACAAACAGATG  
 CTCTACACCAACCAAGTTCCTCAAGAGCAGTGGGTAAAGTTGGAGGACCAGCACGGGAAGCCATATTTCT  
 ACAACCCAGAGGACTCCTCTGTTCAAGTGGGAGTACCCAGGTCCCGATCCCTGCCCCCGAAGTGTCCG  
 CAAATCCAGCCAGGATAGCGACTCCAGCCAGGCTAGCCCTCCAGAGGAGAAGATCAAGACTCTGGAC  
 AAGGCAGGAGTCTCCATCGTACCAAGACGGTAGACAAGGGGAAGAGACTTCGGAAGAAACTGGAGCA  
 CCTCCTGGACAGTGTGGAGGGCGGGTCTGACCTTCTCAAGGACTCAAAGACCTCAGCTGCAGGCGG  
 CCTGAGGCAGCCTTCAAAGCTCTCCACCCTGAATACACAGTCAAGTGAAGGGGCTCTCTCTTTGG  
 GCCCAAGACAAATCCAGCAAGAAGAATGTGTTAGAGCTTCGAAGCCGAGATGGCTCGGAGTACCTGA  
 TCCAGCACGACTCCGAAGCCATCATCAGCACCTGGCACAAGCCATTGCCGAAGGCATCTCAGAGTGTG  
 TGCAGACCTGCTCCAAGGGGAGGAAGGTGAGCCAGCAGTGTGATTTCCGGTCCAGTGAGCGCCTCGGA  
 AGCTGGCGGGAGGAGGATGTACGGCAGAACGCAGCCTCACCTCGCTCAGCCCCGAGGCTGGAAGCG  
 ACTTGAGCAGAGTCCGGCACAAGCTCCGTAATTCCTACAGAGACGACCCACTCTGCAGTCTTTGCGGGA  
 CAAGGGCTACATCAAAGACCAGGTGTTGGATGCGCGCTGGCTCAACTGTGTGAGCGGAGAGGAGCCCT  
 GTGCCACGCTTCGTGCAGCAATGCATCCGACCGTCGAGGCCGGGGCTGGATATTGATGGGCTGTACC  
 GCATCAGTGGGAACCTGGCTACCATCCAGAAGCTGCGCTATAAGGTGGATCACGATGAACGCTTGGACTT  
 GGACGACGGGCGCTGGGAAGACGTCCACGTGATCACGGGTGCCCTGAAGCTCTTCTCCGAGAGCTGCCG  
 GAACCTCTTTCCCTTCTCGCACTTCCACCAGTTCATAGCAGCCATTAAGCTGCAGGACCCAGCCAGC  
 GCAGCCGCTGCGTGCGTGACCTGGTGCGCAGCTGCCAGCTCCCAACCAGGACACGCTCCGACTGCTCAT  
 CCAGCATCTGCGGGTTCATCGAGCATGGCGAGCAGAACCAGTATGACTGTGCAGAAATGTGGCCATAGT  
 TTCGAGCCACGCTGCTGCGGCCAGAGATGGAGGAGGCCAGCATGCCATGACCATGGTGTTCAGAACC  
 AGGTGGTGGAGCTCATCTACATCAGTGCAGGACATCTCCACCGCAC

**ACGGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR224354 representing NM\_133715  
Red=Cloning site Green=Tags(s)


MVDMISKLVRRQSRALRAQVDDPPEPVYANVERQPRATSPRSAAAPPRLSPVWETHDGTGRPYYYNPD  
 TGVTTWESPFETPEGTTSPATSRASVSGESLETEWQYWDDEESRRVFFYNPLTGETAWEDETEELEEDH  
 QEQLQMPSLSRSPGQQRPPTPETDYPELLASYPEEDYSPVGSFSDPGPASPLVAPPWSCQITPDKQM  
 LYTNQFTQEQQWRLEDQHGKPYFYNPEDSSVQWELPQVPIPAPRSVRKSSQSDTPAQASPPEEKIKTLD  
 KAGVLHRTKTVDKGKRLRKKHWSTSWTVLEGGVLTFFKDSKTSAAAGGLRQPSKLSPEYTVELKGASLSW  
 APKDKSSKKNVLELRSRDGSEYLIQHDSEAIISTWHKAI AEGISEL SADLLQGEEGEPSSADFGSSERLG  
 SWREEDVRQNAASPSLSPGGLESLSRVRHKLKFLQRRPTLQSLRDKGYIKDQVFGCALACL CERERS  
 VPRFVQQCIRTVEARGLDIDGLYRISGNLATIQKLRKVDHDERLDLDDGRWEDVHVITGALKLFFRELP  
 EPLFPFSHFHQFIAAIKLDPAQRSRCVRDLVRTLPAPNQDTRLRLLIQHLCRVIEHGEQNRMTVQNVAIV  
 FGPTLLRPEMEEASMPMTMVFNQVVELILHQCADIFPPH

TRTRPLEQKLI SEEDLANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9093\\_e11.zip](https://cdn.origene.com/chromatograms/mm9093_e11.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_133715

**ORF Size:** 2010 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_133715.5](#), [NP\\_598476.2](#)

**RefSeq Size:** 3528 bp

**RefSeq ORF:** 2013 bp

**Locus ID:** 544817


**UniProt ID:** [A2AB59](#)

**Cytogenetics:** 11 E1

**MW:** 75.8 kDa

**Gene Summary:** Rho GTPase-activating protein which may be involved in clathrin-mediated endocytosis. GTPase activators for the Rho-type GTPases act by converting them to an inactive GDP-bound state. Has activity toward CDC42 and RAC1 (By similarity).[UniProtKB/Swiss-Prot Function]

### Product images:



Circular map for MR224354