

Product datasheet for MR224356

Ppm1b (NM_001159498) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ppm1b (NM_001159498) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ppm1b
Synonyms:	PP2CB
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>MR224356 representing NM_001159498 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGTGCATTTTGGATAAACCCAAAACCTGAAAAGCACAATGCTCACGGTGTGGGAATGGTCTGCGTT
ATGGCCTGAGCAGTATGCAAGGATGGAGAGTAGAAATGGAAGATGCACACACAGCTGTTGTGGGTATTCC
TCACGGCTTGGACAACTGGTCGTTTTTTCAGTTTATGACGGTCATGCTGGATCCCGAGTGGCAAATTAC
TGTTCAACACATCTATTAGAACACATCACCACCAATGAAGACTTCAGGGCAGCTGACAAAATCAGGCTCTG
CTCTCGAGCCTTCAGTAGAAAGTGTGAAGACTGGTATCAGGACTGGCTTTTTGAAAATTGATGAATATAT
GCGTAACCTTCTCAGACCTGAGGAACGGGATGGACAGGAGCGGCTCGACTGCAGTGGGCGTGATGGTCTCC
CCTACACACATGTACTTCATCAACTGTGGTGACTCTCGAGCTGTTCTGTGTAGGAACGGACAGGTCTGCT
TTTCTACCCAGGATCACAAACCTTGAATCCAGTGGAGAAGGAGCGCATCCAAAACGCAGGAGGCAGTGT
GATGATCCAGCGTGTGAACGGCTCACTAGCAGTGTCTCGGGCTCTGGGGACTATGACTACAAGTGTGTG
GATGGCAAGGGCCCTACAGAGCAGCTTGTTCCTCAGAGCCTGAGGTTTATGAGATTGTGAGAGCAGAAG
AGGATGAGTTTGTCTTGGCCTGTGATGGGATCTGGGATGTGATGAGCAATGAGGAGCTCTGTGAGTT
TGTTAAGTCTAGGCTTGAGGTGTGCGACGACCTGGAGAATGTGTGCAATTGGGTAGTGGACACTTGTTTA
CATAAGGGAAGTCGAGATAACATGAGTGTGTATTAGTTTGTCTTTCAAATGCCCAAGGTCTCAGAGG
AAGCCGTGAAGAGAGATTGAGATTGGATAAGCACTTGAATCACGGGTTGAAGAAATCATGCAGAAGTC
TGGGGAGGAAGGAATGCCTGATCTTGCCCATGTGATGCGCATTTTGTCTGCAGAAAATATCCCGAATTTA
CCTCCCGGGGAGGTCTCGCTGGCAAGCGCCATGTTATTGAAGCTGTTTATAGTAGACTTAATCCACACA
AAGACAATGATGGGTTTTACCAACCTTCAATTGCATACAGTGACAATGTTTTCTATTA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR224356 representing NM_001159498
 Red=Cloning site Green=Tags(s)


MGAFLDKPKTEKHNAHGAGNGLRYGLSSMQGWRVEMEDAHTAVVGIPHGLDNWSFFAVYDGHAGSRVANY
 CSTHLLLEHITTNEFDRAADKSGSALEPSVESVKTGIRTGFLKIDEYMRNFSDLRNGMDRSGSTAVGVMVS
 PTHMYFINCGDSRAVLCRNGQVCFSTQDHPKPCNPVEKERIQNAGGSVMIQRVNGSLAVSRALGDYDYKCV
 DGKGPTEQLVSPPEVYEVVRAEDEFVVLACDGIWDVMSNEELCEFVKSRLLEVSDDLLENVCNWWVDTCL
 HKGSRDNMSVVLVCFSNAPKVSEEAVKRDELDKHLESRVEEIMQKSGEEGMPDLAHVMRILSAENIPNL
 PPGGGLAGKRHVIEAVYSRLNPHKDNDGFYQPSIAYSNDNVFLL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001159498

ORF Size:	1179 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001159498.2
RefSeq Size:	1580 bp
RefSeq ORF:	1182 bp
Locus ID:	19043
UniProt ID:	P36993
Cytogenetics:	17 55.13 cM
MW:	43.9 kDa
Gene Summary:	Enzyme with a broad specificity. Dephosphorylates PRKAA1 and PRKAA2. Inhibits TBK1-mediated antiviral signaling by dephosphorylating it at 'Ser-172'. Plays an important role in the termination of TNF-alpha-mediated NF-kappa-B activation through dephosphorylating and inactivating IKBKB/IKKB.[UniProtKB/Swiss-Prot Function]