

Product datasheet for **MR224367L4V**

Tm4sf20 (NM_025453) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tm4sf20 (NM_025453) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tm4sf20
Synonyms:	1810018L02Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025453
ORF Size:	678 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224367).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025453.3
RefSeq Size:	1503 bp
RefSeq ORF:	681 bp
Locus ID:	66261
UniProt ID:	Q9CQY8
Cytogenetics:	1 C5



[View online »](#)

Gene Summary:

Polytopic transmembrane protein. Inhibits regulated intramembrane proteolysis (RIP) of CREB3L1, inhibiting its activation and the induction of collagen synthesis. In response to ceramide, which alters TM4SF20 membrane topology, stimulates RIP activation of CREB3L1. Ceramide reverses the direction through which transmembrane helices are translocated into the endoplasmic reticulum membrane during translation of TM4SF20, this mechanism is called 'regulated alternative translocation' (RAT) and regulates the function of the transmembrane protein.[UniProtKB/Swiss-Prot Function]