

Product datasheet for MR224377

Vmn1r219 (NM_134238) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Vmn1r219 (NM_134238) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Vmn1r219
Synonyms:	V1rh13
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR224377 representing NM_134238 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTTTTGCAATTTATTAAGGAAATATTTTTTCTTCATGACTATGGTTGGCATTCTGGGGAACATGT
CTGTTTCTGTGCATTATATGTTCAAGTTGGTGGGAAGCCCTGAGAAGAAACCCATACACCTTATTCTCAC
TCACTTGACTTTTACAAACATCATACTTCTGCTTGCAAAGGATTGCAAAGACAATAACAGTTTTTAGG
TTGAGAACTTCTGGATGACATAGGCTGCAAGATCATTGTTTACCTGGAGAGGGTGGCCCGGGCCTCT
CCATCTGCACCAGCAGTTTCTCACTGTGGTCCAGGCCATCATCATCAGTCCCAGAGCATCTGGGTGGAG
CAGACTCAGGCCAAAGTCTGCATGGCACATCCTTCCATTCTTCTATTCTTTGGATATTCAATGGTTTA
ATAAGTGTGAACCTAATCCATTCGATTAGAAGTATAGGTATGAATACATCACAGCATAAGAATAGTGACA
ACTATTGTTATATGACTAGAAAGTCAGGAAATCAAATGGATTGTTCTCCCTCTCATGGTCTGAGAGA
TGCAGTGTTCAGGGTGCCATGGGAGGGGCCAGTGGCTACATGGTATTTCTTCTCCACAAGCACCACCAG
CAAGTCTCTACTTTCAGAACTCAAAGTCTCTACAGAACTCCCCTGAGCTGAGAGCTGCTCAAAGTG
TCCTCTTCTGATGCTCTGTTTTGTTTTCTTCTACTGGACAGACTGTGCCTTTTCTCTGTTTCTAAGTCT
CTCTTTAGGGGACTATTTGTTGATGGTAAGTATTCATGAATCTGTGACCCTTGGTTATGCAGTTTCCAGC
CCCCTCTGCTGATTCACAGGGATGGACTTTTGACTGAATGTTGCCATTCTCAATGGGAGACATTGAGAA
ACTGTCTCTCTACTATATTTTCAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR224377 representing NM_134238
Red=Cloning site Green=Tags(s)

MVLQFIKEIFFFFMTMVGILGNMSVSVHYMFSSWGSPEKKPIHLILTHLTFTNIILLAKGLQKITVFR
 LRNFLDDIGCKIIVYLERVARGLSICTSSFLTVVQAIISPRASGWSRLRPKSAWHILPFFSFFWIFNGL
 ISVNLHSIRSIGMNTSQHKNSDNYCYMTLESQEIKWIVLPLMVLRDVFGAMGGASGYMVFLHKKHQ
 QVLYFQNSKLLYRTPPELRAAQSVLLMLCFVFFYWTDCAFSLFLSLSLGDYLLMVSIIHESVTLGYAVSS
 PLLLIHRDGLLTCCHSQWETLRNCLSLLYFQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9063_g07.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_134238

ORF Size: 936 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_134238.1](#), [NP_598999.1](#)

RefSeq Size: 939 bp


RefSeq ORF: 939 bp

Locus ID: 171272

Cytogenetics: 13 A3.1

MW: 36 kDa

Product images:



Circular map for MR224377