

Product datasheet for **MR224389L3V**

Slk (NM_001164639) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slk (NM_001164639) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slk
Synonyms:	9A2; AV021402; AW411554; Etk4; mKIAA0204; mSLK; SMAK; Stk2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001164639
ORF Size:	3609 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224389).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001164639.1 , NP_001158111.1
RefSeq Size:	6954 bp
RefSeq ORF:	3609 bp
Locus ID:	20874
UniProt ID:	O54988
Cytogenetics:	19 C3- D1
Gene Summary:	Mediates apoptosis and actin stress fiber dissolution.[UniProtKB/Swiss-Prot Function]



[View online »](#)